



LASER TECHNOLOGY equipment and accessories



THERMO-SPHERIX®

M-RF Thermo Laser Magnetic Rotation Field®



THERMO-INNOVATION

ROTATING MAGNETIC FIELD

10-LASER LIPOLYSIS

BEAUTY TECHNOLOGY WITHOUT SCALPEL

LIFT WITHOUT SURGERY

- ◆ THERMO-FACELIFT
- ◆ THERMO-BODY LIFT

+

LIPOPLASTY

AND BODY THERMO-MODELING

- ◆ SLIMMING ◆ CELLULITE REDUCTION
- ◆ REFIRMING

Rejuvenating volumetric modelling of face, neck and décolleté and body skin anti-aging treatment by thermal rearrangement and reconstruction of elastic fibres, and stimulation of new fibre production, with immediate effects of wrinkle resurfacing, hydrating and skin density improvement.

PREPARATIONS SKIN LIFTING THERMO FASHION



Face treatment effect:

- ◆ Local dynamic face remodelling with non-invasive facelift
- ◆ face contour modelling, chin and local body lift
- ◆ stimulation of elastic fibre production (collagen and elastin) towards restoration of proper natural skin elasticity
- ◆ density reconstruction and skin refirming
- ◆ stimulation of local skin reconstruction in the area of stretch marks, thinning and flabbiness reduction of scars, stretch marks and imperfections,
- ◆ wrinkle resurfacing, also around eyes (reduction of „crow's feet“)
- ◆ reduction of hooded eyes effect
- ◆ cheekbone fat pad reduction

Non-invasive lipoplasty- local body modelling, cellulite reduction and slimming with instantly visible improvement of skin firmness, with effect of correction of skin surface inhomogeneity and improvement of skin quality - colour, moisture.

PREPARATIONS BODY SKIN THERMO FASHION CELLU-ESTHETIC SLIM-ESTHETIC



Body treatment effects:

- ◆ structural body modelling with synchronous refirming and elastic fibre strengthening
- ◆ anti-cellulite effect and neutralizing „orange peel“ skin - special programmes to tackle advanced cellulite problems
- ◆ lipoplasty - reduction of adipose tissue volume and change of its structure
- ◆ skin lift and refirming, density restoration
- ◆ stimulation of metabolic changes in skin and subcutaneous fat tissue - special programmes adjusted to fat tissue thickness
- ◆ smoother and better moisturised skin
- ◆ microcirculation stimulation and detoxification



THERMO-SPHERIX TREATMENTS

THERMO-LIFTING Programme 03



THERMO-MODELLING Programme 07



17

IN-SALON
TREATMENTS

SAFETY
WITHOUT RISK OF SKIN BURN
WITHOUT SCALPEL
PAIN-FREE
NON-INVASIVE
INTUITIVE USER EXPERIENCE

THERMO-SPHERIX®

M-RF Thermo Laser



THERMO-SPHERIX® 2 LASER HEADS + ROTATING MAGNETIC FIELD

- ◆ LASER 1200nm (infrared IR - InfraRed)
- ◆ LASER 940nm (infrared IR - InfraRed)
- ◆ LASER 650nm (red R - Red)
- ◆ Rotating low frequency MAGNETIC FIELD with nominal power 10mT
- ◆ innovative system for a targeted effect of special preparations designed for intensive, non-invasive, thermal skin rejuvenation treatment for men and women to be applied with the innovative equipment THERMO-SPHERIX® M-RF THERMO-LASER, or using another type of cosmetic equipment with heads for which the maximum generated temperature is 70°C

THERMO-SPHERIX
WIDE RANGE
OF FUNCTION:
17 in-salon
treatments

THERMO-SPHERIX® FUNCTION

This is multi-purpose equipment combining the function of two unique devices applied to achieve the following effects: THERMO skin lift rejuvenation and THERMO body modelling. THERMO-SPHERIX® includes as many as 17 special programmes for men and women. The equipment uses the innovative THERMO technology so the treatments are painless, non-invasive and do not carry any risk of skin burn.

THERMO-SPHERIX®:

INNOVATIVE SYNERGY OF MAGNETIC-ELECTRO-PHOTODYNAMIC TECHNOLOGY:

- ◆ patented rotating magnetic field Magnetic Rotation Field® with even controlled temperature distribution in the tissue
- ◆ targeted laser action: 940nm and 1200nm wavelength lasers - the efficient laser energy absorption in adipose cells with effect of thermal accumulation
- ◆ 650nm wavelength laser for temporary decrease of fat cell membrane integrity that enables it to eliminate dissolved adipocytes
- ◆ range of thermo-active preparations Skin Thermo Fashion in applications up to 70 °C
- ◆ intuitive 17- programme control panel with programmed treatment times
- ◆ application in professional beauty and aesthetic medicine salons

Preparations for face thermo-modelling treatment:

SKIN LIFTING THERMO-FASHION KIT Face & Body Skin


SKIN LIFTING THERMO-FASHION Kit for thermo-lifting and thermo-modelling of face skin and local body areas; for application using THERMO-SPHERIX laser equipment and by apparatuses generating a maximum temperature of 70 °C

BB AFTER Resurfacing Day Cream SPF25 UVA/UVB

BB AFTER intensively rejuvenating & regenerative cream with colour balancing pigments, regulating epidermis regrowth, with biomimetic peptide and botanical retinol for skin with visible signs of aging, and damaged, dry skin, recommended also to maintain derm aesthetics treatment effects

AFTER Resurfacing Night Serum Face & Eyes

Night AFTER serum - intensive wrinkle reduction, firming and regeneration of skin of face and around the eyes, regulating epidermis regrowth; contains active anti-wrinkle component Renestril, biomimetic peptide and botanical retinol



Anti-cellulite, slimming and firming preparations for body thermo-modelling treatment:

Tea Body Peel Firming & Shaping Body Skin

Tea peel for slimming and anti-cellulite body care

Slim-Esthetic Body Skin Concentrate

Slim-Esthetic non perfumed slimming gel concentrate

Slim-Esthetic Body Skin Cinnamon Cream

Slim-Esthetic Cinnamon slimming cream

Cellu-Esthetic Body Skin Concentrate

Cellu-Esthetic non perfumed anti cellulite gel concentrate

Cellu-Esthetic Body Skin Grapefruit Cream

Cellu-Esthetic Grapefruit anti-cellulite cream



THERMO-SPHERIX® TECHNICAL SPECIFICATIONS

Power supply: 100-240 V; 50/60 Hz, 20 VA

Lasers R/IR with nominal power: 70 mW

Two laser probes:

1. THERMO-LIFTING for face and local body areas

2. LASER LIPOPLASTY 1200nm+940nm+650nm

Laser safety class: 3R

Rotating magnetic field with nominal power: 10 mT

Areas of operation: labile, contact method, direct probing of surface: LIFTING 4.5 cm², LIPOPLASTY 10.2 m²

Size (width x depth x height): 280mm x 285mm x 25mm

Weight (without accessories) approx. 1.5kg

Warranty: 24 months

www.chantarelle.eu www.chantarelle.pl www.facebook.com/CHANTARELLE.Laboratory.Derm.Aesthetics
Office Poland: biuro@chantarelle.pl tel/fax: +48 22 610 43 10, +48 22 610 44 40 Prof.Cometica, Sierakowska 29, 05-092 Łomianki, Poland



CHANTARELLE
Laboratory Derm Aesthetics

Fundusze Europejskie - dla rozwoju innowacyjnej gospodarki